

# TECARVET 2000

## TECARTHERAPY

### TECHNICAL FEATURES

#### Display

Size 2.6", backlit

#### Power

50W,  $\pm 10\%$

#### Emission Frequency

470KHZ,  $\pm 10\%$

#### Treatment modality

Resistive bipolar

#### Power supply

- Mains
- Battery pack

#### Size


170x220x60 mm

#### Weight

1100 gr

### EQUIPMENT

- 1 Carrying case
- 1 TECARVET 2000 unit
- 1 Bipolar handpiece  $\varnothing$  50 mm
- 1 260 ml gel bottle
- 1 Power supply unit
- 1 Operating manual
- 1 User's guide



# TECARVET 2000

## TECARTHERAPY

Rev. En\_09\_2019

### PROGRAMS LIST

#### DOG / CAT

Osteoarthritis

---

Muscle contracture

---

Tendinopathy

---

Edema